

PCDH19 Antibody



PACO45464

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

Potential calcium-dependent cell-adhesion protein.

Gene ID:

PCDH19

Uniprot

Q8TAB3

Synonyms:

Protocadherin-19, PCDH19, KIAA1313

Immunogen:

Recombinant Human Protocadherin-19 protein (568-678AA).


Storage:

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Immunohistochemistry of paraffin-embedded human brain tissue using PACO45464 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human glioma using PACO45464 at dilution of 1:100.