

TGFBR3 Antibody



PACO45488

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

Binds to TGF-beta. Could be involved in capturing and retaining TGF-beta for presentation to the signaling receptors.

Gene ID:

TGFBR3

Uniprot

Q03167

Synonyms:

Transforming growth factor beta receptor type 3 (TGF-beta receptor type 3) (TGF-3) (Betaglycan) (Transforming growth factor beta receptor III) (TGF-beta receptor type III), TGFBR3


Immunogen:

Recombinant Human Transforming growth factor β receptor type 3 protein (140-340AA).

Storage:

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Immunohistochemistry of paraffin-embedded human kidney tissue using PACO45488 at dilution of 1:100.