

GBF1 Antibody



PACO45526

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC, IF

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200,
IF:1:50-1:200

Protein Background:**Gene ID:**

GBF1

Uniprot

Q92538

Synonyms:

Golgi-specific brefeldin A-resistance guanine nucleotide exchange factor 1 (BFA-resistant GEF 1), GBF1, KIAA0248


Immunogen:

Recombinant Human Golgi-specific brefeldin A-resistance guanine nucleotide exchange factor 1 protein (1-85AA).


Storage:

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.


Product Images



Immunohistochemistry of paraffin-embedded human kidney tissue using PACO45526 at dilution of 1:100.



Immunofluorescent analysis of HeLa cells using PACO45526 at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L).



Immunohistochemistry of paraffin-embedded human gastric cancer using PACO45526 at dilution of 1:100.