

USP2 Antibody



PACO45631

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:**Gene ID:**

USP2

Uniprot

O75604

Synonyms:

Ubiquitin carboxyl-terminal hydrolase 2 (EC 3.4.19.12) (41 kDa ubiquitin-specific protease) (Deubiquitinating enzyme 2) (Ubiquitin thioesterase 2) (Ubiquitin-specific-processing protease 2), USP2, UBP41


Immunogen:

Recombinant Human Ubiquitin carboxyl-terminal hydrolase 2 protein (200-430AA).


Storage:

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Immunohistochemistry of paraffin-embedded human brain tissue using PACO45631 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human glioma using PACO45631 at dilution of 1:100.