

USP28 Antibody



PACO45632

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

Gene ID:

USP28

Uniprot

Q96RU2

Synonyms:

Ubiquitin carboxyl-terminal hydrolase 28 (EC 3.4.19.12) (Deubiquitinating enzyme 28) (Ubiquitin thioesterase 28) (Ubiquitin-specific-processing protease 28), USP28, KIAA1515


Immunogen:

Recombinant Human Ubiquitin carboxyl-terminal hydrolase 28 protein (430-580AA).

Storage:

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Immunohistochemistry of paraffin-embedded human ovarian cancer using PACO45632 at dilution of 1:100.