

TMEFF2 Antibody



PACO45643

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

May be a survival factor for hippocampal and mesencephalic neurons. The shedded form up-regulates cancer cell proliferation, probably by promoting ERK1/2 phosphorylation.

Gene ID:

TMEFF2

Uniprot

Q9UIK5

Synonyms:

Tomoregulin-2 (TR-2) (Hyperplastic polyposis protein 1) (Transmembrane protein with EGF-like and two follistatin-like domains), TMEFF2, HPP1 TENB2 TPEF


Immunogen:

Recombinant Human Tomoregulin-2 protein (41-320AA).

Storage:

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Immunohistochemistry of paraffin-embedded human melanoma using PACO45643 at dilution of 1:100.