

APOBEC3G Antibody



PACO45670

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC, IF

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200,
IF:1:50-1:200

Protein Background:

Gene ID:

APOBEC3G

Uniprot

Q9HC16

Synonyms:

DNA dC->dU-editing enzyme APOBEC-3G (EC 3.5.4. -) (APOBEC-related cytidine deaminase) (APOBEC-related protein) (ARCD) (APOBEC-related protein 9) (ARP-9) (CEM-15) (CEM15) (Deoxycytidine deaminase) (A3G), APOBEC3G


Immunogen:

Recombinant Human DNA dC->dU-editing enzyme APOBEC-3G protein (1-384AA).


Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Immunohistochemistry of paraffin-embedded human liver tissue using PACO45670 at dilution of 1:100.



Immunofluorescent analysis of MCF-7 cells using PACO45670 at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L).