

APOBEC3G Antibody, Biotin conjugated



PACO45673

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

APOBEC3G

Uniprot

Q9HC16

Synonyms:

DNA dC->dU-editing enzyme APOBEC-3G (EC 3.5.4. -) (APOBEC-related cytidine deaminase) (APOBEC-related protein) (ARCD) (APOBEC-related protein 9) (ARP-9) (CEM-15) (CEM15) (Deoxycytidine deaminase) (A3G), APOBEC3G

Immunogen:

Recombinant Human DNA dC->dU-editing enzyme APOBEC-3G protein (1-384AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A