

# BUB1 Antibody



PACO45698

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, IHC

**Recommended dilutions:**

ELISA:1:2000-1:10000, IHC:1:20-1:200

**Protein Background:****Gene ID:**

BUB1

**Uniprot**

O43683

**Synonyms:**

Mitotic checkpoint serine/threonine-protein kinase BUB1 (hBUB1) (EC 2.7.11.1) (BUB1A), BUB1, BUB1L

**Immunogen:**


Recombinant Human Mitotic checkpoint serine/threonine-protein kinase BUB1 protein (256-577AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---



Immunohistochemistry of paraffin-embedded human gastric cancer using PACO45698 at dilution of 1:100.