

CDSN Antibody, FITC conjugated



PACO45732

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Important for the epidermal barrier integrity.

Gene ID:

CDSN

Uniprot

Q15517

Synonyms:

Corneodesmosin (S protein), CDSN

Immunogen:

Recombinant Human Corneodesmosin protein (40-229AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A