

CEP290 Antibody, Biotin conjugated



PACO45741

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:****Gene ID:**

CEP290

Uniprot

O15078

Synonyms:

Centrosomal protein of 290 kDa (Cep290) (Bardet-Biedl syndrome 14 protein) (Cancer/testis antigen 87) (CT87) (Nephrocystin-6) (Tumor antigen se2-2), CEP290, BBS14 KIAA0373 NPHP6

Immunogen:

Recombinant Human Centrosomal protein of 290 kDa protein (86-156AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A