

CES3 Antibody, FITC conjugated



PACO45744

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Involved in the detoxification of xenobiotics and in the activation of ester and amide prodrugs. Shows low catalytic efficiency for hydrolysis of CPT-11 (7-ethyl-10-[4-(1-piperidino)-1-piperidino]-carbonyloxycamptothecin), a prodrug for camptothecin used in cancer therapeutics.

Gene ID:

CES3

Uniprot

Q6UWW8

Synonyms:

Carboxylesterase 3 (EC 3.1.1.1) (Liver carboxylesterase 31 homolog), CES3

Immunogen:

Recombinant Human Carboxylesterase 3 protein (65-233AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A