

DNMT1 Antibody



PACO45774

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC, IF

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200,
IF:1:50-1:200

Protein Background:**Gene ID:**

DNMT1

Uniprot

P26358

Synonyms:

DNA (cytosine-5)-methyltransferase 1 (Dnmt1) (EC 2.1.1.37) (CXXC-type zinc finger protein 9) (DNA methyltransferase Hsal) (DNA MTase Hsal) (M. Hsal) (MCMT), DNMT1, AIM CXXC9 DNMT


Immunogen:

Recombinant Human DNA (cytosine-5)-methyltransferase 1 protein (1288-1632AA).


Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Immunohistochemistry of paraffin-embedded human colon cancer using PACO45774 at dilution of 1:100.



Immunofluorescent analysis of HepG2 cells using PACO45774 at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L).