

# ELK1 Antibody, Biotin conjugated



PACO45801

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

Transcription factor that binds to purine-rich DNA sequences. Forms a ternary complex with SRF and the ETS and SRF motifs of the serum response element (SRE) on the promoter region of immediate early genes such as FOS and IER2. Induces target gene transcription upon JNK-signaling pathway stimulation.

**Gene ID:**

ELK1

**Uniprot**

P19419

**Synonyms:**

ETS domain-containing protein Elk-1, ELK1

**Immunogen:**

Recombinant Human ETS domain-containing protein Elk-1 protein (247-403AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A