

# EPM2A Antibody, HRP conjugated



PACO45807

---

## Product Information

<b>Size:</b>	<b>Protein Background:</b>
50ug	negative regulation of TOR signaling, positive regulation of macroautophagy.
<b>Reactivity:</b>	<b>Gene ID:</b>
Human	EPM2A
<b>Source:</b>	<b>Uniprot</b>
Rabbit	B3EWF7
<b>Isotype:</b>	<b>Synonyms:</b>
IgG	Laforin, isoform 9, EPM2A
<b>Applications:</b>	<b>Immunogen:</b>
ELISA	Recombinant Human Laforin, isoform 9 protein (1-331AA).
<b>Recommended dilutions:</b>	<b>Storage:</b>
	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A