

EPM2A Antibody, FITC conjugated



PACO45808

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

negative regulation of TOR signaling, positive regulation of macroautophagy.

Gene ID:

EPM2A

Uniprot

B3EWF7

Synonyms:

Laforin, isoform 9, EPM2A

Immunogen:

Recombinant Human Laforin, isoform 9 protein (1-331AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A