

EPM2A Antibody, Biotin conjugated



PACO45809

Product Information

Size:	Protein Background:
50ug	negative regulation of TOR signaling, positive regulation of macroautophagy.
Reactivity:	Gene ID:
Human	EPM2A
Source:	Uniprot
Rabbit	B3EWF7
Isotype:	Synonyms:
IgG	Laforin, isoform 9, EPM2A
Applications:	Immunogen:
ELISA	Recombinant Human Laforin, isoform 9 protein (1-331AA).
Recommended dilutions:	Storage:
	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A