

FOXG1 Antibody, HRP conjugated



PACO45839

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Transcription repression factor which plays an important role in the establishment of the regional subdivision of the developing brain and in the development of the telencephalon.

Gene ID:

FOXG1

Uniprot

P55316

Synonyms:

Forkhead box protein G1 (Brain factor 1) (BF-1) (BF1) (Brain factor 2) (BF-2) (BF2) (hBF-2) (Forkhead box protein G1A) (Forkhead box protein G1B) (Forkhead box protein G1C) (Forkhead-related protein FKHL1) (HFK1) (Forkhead-related protein FKHL2) (HFK2) (Forkhead-related protein FKHL3) (HFK3), FOXG1, FKH2 FKHL1 FKHL2 FKHL3 FKHL4 FOXG1A FOXG1B FOXG1C

Immunogen:

Recombinant Human Forkhead box protein G1 protein (183-292AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A