

INPP4B Antibody, Biotin conjugated



PACO45937

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Catalyzes the hydrolysis of the 4-position phosphate of phosphatidylinositol 3,4-bisphosphate, inositol 1,3,4-trisphosphate and inositol 1,4-bisphosphate.

Gene ID:

INPP4B

Uniprot

O15327

Synonyms:

Type II inositol 3,4-bisphosphate 4-phosphatase (EC 3.1.3.66) (Inositol polyphosphate 4-phosphatase type II), INPP4B

Immunogen:

Recombinant Human Type II inositol 3,4-bisphosphate 4-phosphatase protein (1-53AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A