

IRF5 Antibody, FITC conjugated



PACO45948

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Transcription factor involved in the induction of interferons IFNA and INFB and inflammatory cytokines upon virus infection. Activated by TLR7 or TLR8 signaling.

Gene ID:

IRF5

Uniprot

Q13568

Synonyms:

Interferon regulatory factor 5 (IRF-5), IRF5

Immunogen:

Recombinant Human Interferon regulatory factor 5 protein (214-395AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A