

PACO45974

Product Information

Size:

50ug

Reactivity:

Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IF

Recommended dilutions:

ELISA:1:2000-1:10000, IF:1:50-1:200

Protein Background:

Galectin that binds lactose and a related range of sugars.

Gene ID:

Lgals4

Uniprot

P38552

Synonyms:

Galectin-4 (Gal-4) (L-36 lactose-binding protein) (L36LBP) (Lactose-binding lectin 4), Lgals4


Immunogen:

Recombinant Rat Galectin-4 protein (1-324AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Immunofluorescent analysis of HepG2 cells using PACO45974 at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L).