

LRRIQ3 Antibody, FITC conjugated



PACO45992

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:****Gene ID:**

LRRIQ3

Uniprot

A6PVS8

Synonyms:

Leucine-rich repeat and IQ domain-containing protein 3 (Leucine-rich repeat-containing protein 44), LRRIQ3, LRRC44

Immunogen:

Recombinant Human Leucine-rich repeat and IQ domain-containing protein 3 protein (538-607AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A