

MPL Antibody, Biotin conjugated



PACO46037

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Receptor for thrombopoietin. May represent a regulatory molecule specific for TPO-R-dependent immune responses.

Gene ID:

MPL

Uniprot

P40238

Synonyms:

Thrombopoietin receptor (TPO-R) (Myeloproliferative leukemia protein) (Proto-oncogene c-Mpl) (CD antigen CD110), MPL, TPOR

Immunogen:

Recombinant Human Thrombopoietin receptor protein (26-491AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A