

MYH3 Antibody, HRP conjugated



PACO46043

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Muscle contraction.

Gene ID:

MYH3

Uniprot

P11055

Synonyms:

Myosin-3 (Muscle embryonic myosin heavy chain) (Myosin heavy chain 3) (Myosin heavy chain, fast skeletal muscle, embryonic) (SMHCE), MYH3

Immunogen:

Recombinant Human Myosin-3 protein (1-217AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A