

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC, IF

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200,
IF:1:50-1:200

Protein Background:

High affinity receptor for semaphorins 3C, 3F, VEGF-165 and VEGF-145 isoforms of VEGF, and the PLGF-2 isoform of PGF.

Gene ID:

NRP2

Uniprot

O60462

Synonyms:

Neuropilin-2 (Vascular endothelial cell growth factor 165 receptor 2), NRP2, VEGF165R2


Immunogen:

Recombinant Human Neuropilin-2 protein (60-209AA).


Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4


Product Images



Immunohistochemistry of paraffin-embedded human colon cancer using PACO46090 at dilution of 1:100.



Immunofluorescent analysis of HepG2 cells using PACO46090 at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L).



Immunohistochemistry of paraffin-embedded human skeletal muscle tissue using PACO46090 at dilution of 1:100.