

PCF11 Antibody, Biotin conjugated



PACO46101

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Component of pre-mRNA cleavage complex II.

Gene ID:

PCF11

Uniprot

O94913

Synonyms:

Pre-mRNA cleavage complex 2 protein Pcf11 (Pre-mRNA cleavage complex II protein Pcf11), PCF11, KIAA0824

Immunogen:

Recombinant Human Pre-mRNA cleavage complex 2 protein Pcf11 protein (409-604AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A