

PACO46118

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, IHC, IF

**Recommended dilutions:**

ELISA:1:2000-1:10000, IHC:1:20-1:200,  
IF:1:50-1:200

**Protein Background:**

**Gene ID:**

FAU

**Uniprot**

P62861

**Synonyms:**

40S ribosomal protein S30 (Small ribosomal subunit protein eS30), FAU

**Immunogen:**


Recombinant Human 40S ribosomal protein S30 protein (1-59AA).

**Storage:**


Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images


---



Immunohistochemistry of paraffin-embedded human endometrial cancer using PACO46118 at dilution of 1:100.



Immunofluorescent analysis of HeLa cells using PACO46118 at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L).



Immunohistochemistry of paraffin-embedded human colon cancer using PACO46118 at dilution of 1:100.