

FAU Antibody, FITC conjugated



PACO46120

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:****Gene ID:**

FAU

Uniprot

P62861

Synonyms:

40S ribosomal protein S30 (Small ribosomal subunit protein eS30), FAU

Immunogen:

Recombinant Human 40S ribosomal protein S30 protein (1-59AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A