

PRSS12 Antibody, Biotin conjugated



PACO46145

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Plays a role in neuronal plasticity and the proteolytic action may subserve structural reorganizations associated with learning and memory operations.

Gene ID:

PRSS12

Uniprot

P56730

Synonyms:

Neurotrypsin (EC 3.4.21) (Leydin) (Motopsin) (Serine protease 12), PRSS12

Immunogen:

Recombinant Human Neurotrypsin protein (631-874AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A