

# RHD Antibody, FITC conjugated



PACO46168

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

May be part of an oligomeric complex which is likely to have a transport or channel function in the erythrocyte membrane.

**Gene ID:**

RHD

**Uniprot**

Q02161

**Synonyms:**

Blood group Rh(D) polypeptide (RHXIII) (Rh polypeptide 2) (RhPII) (Rhesus D antigen) (CD antigen CD240D), RHD

**Immunogen:**

Recombinant Human Blood group Rh(D) polypeptide protein (32-76AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A