

SCG2 Antibody, HRP conjugated



PACO46195

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Secretogranin-2 is a neuroendocrine secretory granule protein, which is the precursor for biologically active peptides.

Gene ID:

SCG2

Uniprot

P13521

Synonyms:

Secretogranin-2 (Chromogranin-C) (Secretogranin II) (SgII) [Cleaved into: Secretoneurin (SN); Manserin], SCG2, CHGC

Immunogen:

Recombinant Human Secretogranin-2 protein (31-617AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A