

SOX14 Antibody, HRP conjugated



PACO46223

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Acts as a negative regulator of transcription.

Gene ID:

SOX14

Uniprot

O95416

Synonyms:

Transcription factor SOX-14 (Protein SOX-28), SOX14, SOX28

Immunogen:

Recombinant Human Transcription factor SOX-14 protein (1-137AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A