

# TAF6 Antibody, HRP conjugated



PACO46263

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

TAFs are components of the transcription factor IID (TFIID) complex, PCAF histone acetylase complex and TBP-free TAFII complex (TFTC). TFIID is multimeric protein complex that plays a central role in mediating promoter responses to various activators and repressors.

**Gene ID:**

TAF6

**Uniprot**

P49848

**Synonyms:**

Transcription initiation factor TFIID subunit 6 (RNA polymerase II TBP-associated factor subunit E) (Transcription initiation factor TFIID 70 kDa subunit) (TAF(II)70) (TAFII-70) (TAFII70) (Transcription initiation factor TFIID 80 kDa subunit) (TAF(II)80) (TAFII-80) (TAFII80), TAF6, TAF2E TAFII70

**Immunogen:**

Recombinant Human Transcription initiation factor TFIID subunit 6 protein (627-677AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A