

Tph1 Antibody



PACO46294

Product Information

Size:

50ug

Reactivity:

Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IF

Recommended dilutions:

ELISA:1:2000-1:10000, IF:1:50-1:200

Protein Background:

cytosol, indolalkylamine biosynthetic process.

Gene ID:

Tph1

Uniprot

P17532

Synonyms:

Tryptophan 5-hydroxylase 1 (EC 1.14.16.4) (Tryptophan 5-monoxygenase 1), Tph1, Tph


Immunogen:

Recombinant Mouse Tryptophan 5-hydroxylase 1 protein (1-447AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Immunofluorescent analysis of HepG2 cells using PACO46294 at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L).