

# ZNF408 Antibody, FITC conjugated



PACO46340

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

May be involved in transcriptional regulation.

**Gene ID:**

ZNF408

**Uniprot**

Q9H9D4

**Synonyms:**

Zinc finger protein 408 (PR domain zinc finger protein 17), ZNF408, PFM14 PRDM17

**Immunogen:**

Recombinant Human Zinc finger protein 408 protein (317-635AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A