

ZNF408 Antibody, Biotin conjugated



PACO46341

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

May be involved in transcriptional regulation.

Gene ID:

ZNF408

Uniprot

Q9H9D4

Synonyms:

Zinc finger protein 408 (PR domain zinc finger protein 17), ZNF408, PFM14 PRDM17

Immunogen:

Recombinant Human Zinc finger protein 408 protein (317-635AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A