

ANGPTL1 Antibody



PACO46382

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

extracellular exosome, extracellular space, receptor binding.

Gene ID:

ANGPTL1

Uniprot

O95841

Synonyms:

Angiotensin-related protein 1 (Angiotensin-3) (ANG-3) (Angiotensin-like protein 1), ANGPTL1, ANG3 ANGPT3 ARP1


Immunogen:

Recombinant Human Angiotensin-related protein 1 protein (24-289AA).


Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Immunohistochemistry of paraffin-embedded human kidney tissue using PACO46382 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human testis tissue using PACO46382 at dilution of 1:100.