

ANGPTL1 Antibody, HRP conjugated



PACO46383

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

extracellular exosome, extracellular space, receptor binding.

Gene ID:

ANGPTL1

Uniprot

O95841

Synonyms:

Angiotensin-converting enzyme 2 (ACE2) (Angiotensin-converting enzyme 2-like protein 1), ANGPTL1, ANG3 ANGPT3 ARP1

Immunogen:

Recombinant Human Angiotensin-converting enzyme 2 protein (24-289AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A