

# ANGPTL1 Antibody, FITC conjugated



PACO46384

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

extracellular exosome, extracellular space, receptor binding.

**Gene ID:**

ANGPTL1

**Uniprot**

O95841

**Synonyms:**

Angiopoietin-related protein 1 (Angiopoietin-3) (ANG-3) (Angiopoietin-like protein 1), ANGPTL1, ANG3 ANGPT3 ARP1

**Immunogen:**

Recombinant Human Angiopoietin-related protein 1 protein (24-289AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A