

# MLH1 Antibody, HRP conjugated



PACO46411

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

MLH1

**Uniprot**

P40692

**Synonyms:**

DNA mismatch repair protein Mlh1 (MutL protein homolog 1), MLH1, COCA2

**Immunogen:**

Recombinant Human DNA mismatch repair protein Mlh1 protein (350-519AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A