

PACO46422

---

## Product Information

**Size:**

50ug

**Reactivity:**

Arabidopsis thaliana

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

Produces ATP from ADP in the presence of a proton gradient across the membrane. The catalytic sites are hosted primarily by the beta subunits.

**Gene ID:**

atpB

**Uniprot**

P19366

**Synonyms:**

ATP synthase subunit beta, chloroplastic (EC 3.6.3.14) (ATP synthase F1 sector subunit beta) (F-ATPase subunit beta), atpB

**Immunogen:**

Recombinant Arabidopsis thaliana ATP synthase subunit beta, chloroplastic protein (92-393AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A