

PSBO Antibody, HRP conjugated



PACO46435

Product Information

Size:

50ug

Reactivity:

Spinacia oleracea

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Stabilizes the manganese cluster which is the primary site of water splitting. Binds GTP after preillumination of photosystem II core complex. This binding is inhibited by DCMU.

Gene ID:

PSBO

Uniprot

P12359

Synonyms:

Oxygen-evolving enhancer protein 1, chloroplastic (OEE1) (33 kDa subunit of oxygen evolving system of photosystem II) (33 kDa thylakoid membrane protein) (OEC 33 kDa subunit), PSBO

Immunogen:

Recombinant Spinacia oleracea Oxygen-evolving enhancer protein 1, chloroplastic protein (63-232AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A