

ISPD Antibody, Biotin conjugated



PACO46457

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Required for protein O-linked mannosylation. Probably acts as a nucleotidyltransferase involved in synthesis of a nucleotide sugar. Required for dystroglycan O-mannosylation.

Gene ID:

ISPD

Uniprot

A4D126

Synonyms:

D-ribitol-5-phosphate cytidyltransferase (EC 2.7.7.40) (2-C-methyl-D-erythritol 4-phosphate cytidyltransferase-like protein) (Isoprenoid synthase domain-containing protein) (hISPD), ISPD

Immunogen:

Recombinant Human Isoprenoid synthase domain-containing protein (1-264AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A