

ARHGAP42 Antibody, HRP conjugated



PACO46459

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

May influence blood pressure by functioning as a GTPase-activating protein for RHOA in vascular smooth muscle.

Gene ID:

ARHGAP42

Uniprot

A6NI28

Synonyms:

Rho GTPase-activating protein 42 (Rho GTPase-activating protein 10-like) (Rho-type GTPase-activating protein 42), ARHGAP42, GRAF3

Immunogen:

Recombinant Human Rho GTPase-activating protein 42 protein (541-728AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A