

# SMIM1 Antibody, FITC conjugated



PACO46464

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

Regulator of red blood cells formation.

**Gene ID:**

SMIM1

**Uniprot**

B2RUZ4

**Synonyms:**

Small integral membrane protein 1 (Vel blood group antigen), SMIM1

**Immunogen:**

Recombinant Human Small integral membrane protein 1 protein (1-46AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A