

FSTL1 Antibody, FITC conjugated



PACO46480

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

May modulate the action of some growth factors on cell proliferation and differentiation. Binds heparin.

Gene ID:

FSTL1

Uniprot

Q12841

Synonyms:

Follistatin-related protein 1 (Follistatin-like protein 1), FSTL1, FRP

Immunogen:

Recombinant Human Follistatin-related protein 1 protein (176-285AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A