

# FSTL1 Antibody, Biotin conjugated



PACO46481

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

May modulate the action of some growth factors on cell proliferation and differentiation. Binds heparin.

**Gene ID:**

FSTL1

**Uniprot**

Q12841

**Synonyms:**

Follistatin-related protein 1 (Follistatin-like protein 1), FSTL1, FRP

**Immunogen:**

Recombinant Human Follistatin-related protein 1 protein (176-285AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A