

VGLL4 Antibody, HRP conjugated



PACO46487

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

May act as a specific coactivator for the mammalian TEFs.

Gene ID:

VGLL4

Uniprot

Q14135

Synonyms:

Transcription cofactor vestigial-like protein 4 (Vgl-4), VGLL4, KIAA0121

Immunogen:

Recombinant Human Transcription cofactor vestigial-like protein 4 protein (43-277AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A