

ZRSR2P1 Antibody, HRP conjugated



PACO46499

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:****Gene ID:**

ZRSR2P1

Uniprot

Q15695

Synonyms:

Putative U2 small nuclear ribonucleoprotein auxiliary factor 35 kDa subunit-related protein 1 (CCCH type zinc finger, RNA-binding motif and serine/arginine rich protein 1) (U2(RNU2) small nuclear RNA auxiliary factor 1-like 1), ZRSR2P1, U2AF1-RS1, U2AF1L1, U2AF1P, U2AF1RS1, U2AFBPL, ZRSR1

Immunogen:

Recombinant Human Putative U2 small nuclear ribonucleoprotein auxiliary factor 35 kDa subunit-related protein 1 protein (197-372AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A