

LPAL2 Antibody, FITC conjugated



PACO46520

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:****Gene ID:**

LPAL2

Uniprot

Q16609

Synonyms:

Putative apolipoprotein(a)-like protein 2 (Apo(a)-like protein 2) (Lp(a)-liker protein 2) (Apolipoprotein a-related gene C protein) (Apo(a)rg-C), LPAL2, APOARGC

Immunogen:

Recombinant Human Putative apolipoprotein(a)-like protein 2 protein (52-124AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A